



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

HEY2 Antibody [#C11818]

Catalog Number: C11818

Concentration: 1mg/ml

Swiss-Prot No.: Q9UBP5

Other Names: Hairy/enhancer-of-split related with YRPW motif protein 2, Cardiovascular helix-loop-helix factor 1, hCHF1, Hairy and enhancer of split-related protein 2, HESR-2, Hairy-related transcription factor 2, hHRT2, HES-related repressor protein 2, Protein gridlock homolog, HEY2, CHF1, GRL, HERP, HERP1, HRT2

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against

synthesized peptide derived from internal of human HEY2.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: HEY2 antibody detects endogenous levels of total HEY2 protein.

Reactivity: Human, Mouse

Applications: WB: 1:500~1:3000 ELISA: 1:1000

References:

Chin M.T., J. Biol. Chem. 275:6381-6387(2000).

Nakagawa O., Proc. Natl. Acad. Sci. U.S.A. 97:13655-13660(2000).

Iso T., Mol. Cell. Biol. 21:6071-6079(2001).

