



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

ZC3H7B Antibody [#C11837]

Catalog Number: C11837

Concentration: 1mg/ml

Swiss-Prot No.: Q9UGR2

Other Names: KIAA1031, ROXN, Rotavirus 'X' associated non-structural protein, Z3H7B, ZC3H7B, zinc finger CCCH-type containing 7B

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human

ZC3H7B.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: ZC3H7B antibody detects endogenous levels of total ZC3H7B protein.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:10000

References:

Ota T., Nat. Genet. 36:40-45(2004).

Rush J., Nat. Biotechnol. 23:94-101(2005).

Kikuno R., DNA Res. 6:197-205(1999).

