



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

PLMN (heavy chain A short form, Cleaved-Val98) Antibody [#L0370]

Catalog Number: L0370

Concentration: 1mg/ml

Swiss-Prot No.: P00747

Other Names: PLG; PLMN

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from human PLMN.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: PLMN (heavy chain A short form,

Cleaved-Val98) antibody detects endogenous levels of fragment of activated PLMN resulting from cleavage adjacent to Val98.

Reactivity: Human, Mouse, Rat


Applications: WB: 1:500~1:3000 ELISA: 1:1000

References:

Muller D., Bioche Petersen T.E., J. Biol. Chem. 265:6104-6111(1990).

Forsgren M., FEBS Lett. 213:254-260(1987).

Browne M.J., Submitted (OCT-1991) to the EMBL/GenBank/DDBJ databases.



For Research Use Only