



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

CNOT7 Antibody [#C11858]

Catalog Number: C11858

Concentration: 1mg/ml

Swiss-Prot No.: Q9UIV1

Other Names: BTG1 binding factor 1, BTG1-binding factor 1, CCR4-NOT transcription complex, 7, CCR4-associated factor 1, carbon catabolite repressor protein (CCR4)-associative factor 1

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from N-terminal of human

CNOT7.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: CNOT7 antibody detects endogenous levels of total CNOT7 protein.

Reactivity: Human, Mouse

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Bogdan J.A., Biochem. J. 336:471-481(1998).

The MGC Project Team; Genome Res. 14:2121-2127(2004).

