



**Web:** [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
**Order:** [order@dijibio.com](mailto:order@dijibio.com)  
**Tel:** (86) 519 8805 0026

## STAG3 Antibody [#C11859]

**Catalog Number:** C11859

**Concentration:** 1mg/ml

**Swiss-Prot No.:** Q9UJ98

**Other Names:** Cohesin subunit SA-3, Stromal antigen 3, Stromalin 3, SCC3 homolog 3, STAG3

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from C-terminal of human STAG3.

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** STAG3 antibody detects endogenous levels of total STAG3 protein.

**Reactivity:** Human

**Applications:** WB: 1:500~1:3000 ELISA: 1:10000

**References:**

Pezzi N., FASEB J. 14:581-592(2000).

The German cDNA consortium; Submitted (JUL-2002) to the EMBL/GenBank/DDBJ databases.

Prieto I., EMBO Rep. 3:543-550(2002).

