



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

ZNF225 Antibody [#C11865]

Catalog Number: C11865

Concentration: 1mg/ml

Swiss-Prot No.: Q9UK10

Other Names: Zinc finger protein 225, ZNF225, ZN225

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from C-terminal of human ZNF225.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: ZNF225 antibody detects endogenous levels of total ZNF225 protein.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Shannon M., Genome Res. 13:1097-1110(2003).

ood J., Nature 428:529-535(2004).

The MGC Project Team; Genome Res. 14:2121-2127(2004).

