



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

IKZF3 Antibody [#C11873]

Catalog Number: C11873

Concentration: 1mg/ml

Swiss-Prot No.: Q9UKT9

Other Names: AIOL, IKZF3, ZNFN1A3, zinc finger protein Aiolos

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human IKZF3.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: IKZF3 antibody detects endogenous levels of total IKZF3 protein.

Reactivity: Human, Mouse

Applications: WB: 1:500~1:3000 ELISA: 1:10000

References:

Hosokawa Y., Genomics 61:326-329(1999).

Liippo J.P., Eur. J. Immunol. 31:3469-3474(2001).

Perdomo J., J. Biol. Chem. 275:38347-38354(2000).

