



**Web:** [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
**Order:** [order@dijibio.com](mailto:order@dijibio.com)  
**Tel:** (86) 519 8805 0026

## ELF5 Antibody [#C11875]

**Catalog Number:** C11875

**Concentration:** 1mg/ml

**Swiss-Prot No.:** Q9UKW6

**Other Names:** ETS-related transcription factor Elf-5, E74-like factor 5, Epithelium-specific Ets transcription factor 2, ESE-2, Epithelium-restricted ESE-1-related Ets factor, ELF5, ESE2

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against

synthesized peptide derived from internal of human ELF5.

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** ELF5 antibody detects endogenous levels of total ELF5 protein.

**Reactivity:** Human, Mouse

**Applications:** WB: 1:500~1:3000 ELISA: 1:20000

**References:**

Taylor T.D., Nature 440:497-500(2006).

Ota T., Nat. Genet. 36:40-45(2004).

Oettgen P., J. Biol. Chem. 274:29439-29452(1999).

