



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

Mucin-14 Antibody [#C11879]

Catalog Number: C11879

Concentration: 1mg/ml

Swiss-Prot No.: Q9ULC0

Other Names: Endomucin [Precursor], Endomucin-2, Mucin-14, Gastric cancer antigen Ga34, EMCN, EMCN2, MUC14

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human Mucin-14.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: Mucin-14 antibody detects endogenous levels of total Mucin-14 protein.

Reactivity: Human (Identities = 100%, Positives = 100%);
Mouse (Identities = 100%, Positives = 100%);
Rat (Identities = 100%, Positives = 100%)

Applications: IF: 1:100~1:500 ELISA: 1:1000

References:

Kinoshita M., FEBS Lett. 499:121-126(2001).

Kuhn A., J. Invest. Dermatol. 119:1388-1393(2002).

The MGC Project Team, Genome Res. 14:2121-2127(2004).

