



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

PKCB1 Antibody [#C11882]

Catalog Number: C11882

Concentration: 1mg/ml

Swiss-Prot No.: Q9ULU4

Other Names: CTCL tumor antigen se14-3, Cutaneous T-cell lymphoma associated antigen se14-3, KIAA1125, PKCB, PKCB1, PRKCBP1, Protein kinase C binding protein 1, ZMYND8, Zinc finger MYND domain containing protein 8

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human

PKCB1.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: PKCB1 antibody detects endogenous levels of total PKCB1 protein.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:1000

References:

Fossey S.C., Mamm. Genome 11:919-925(2000).

Eichmueller S., Proc. Natl. Acad. Sci. U.S.A. 98:629-634(2001).

Deloukas P., Nature 414:865-871(2001).

