



**Web:** [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
**Order:** [order@dijibio.com](mailto:order@dijibio.com)  
**Tel:** (86) 519 8805 0026

## CFAB Bb (Cleaved-Lys260) Antibody [#L0234]

**Catalog Number:** L0234

**Concentration:** 1mg/ml

**Swiss-Prot No.:** P00751

**Other Names:** Complement factor B [Precursor], EC 3.4.21.47, C3/C5 convertase, Properdin factor B, Glycine-rich beta glycoprotein, GBG, PBF2, Complement factor B Ba fragment, Complement factor B Bb fragment, CFB, BF, BFD

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against

synthesized peptide derived from human CFAB Bb.

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** CFAB Bb (Cleaved-Lys260) antibody detects endogenous levels of fragment of activated CFAB Bb resulting from cleavage adjacent to Lys260.

**Reactivity:** Human, Mouse

**Applications:** WB: 1:500~1:3000 ELISA: 1:20000

**References:**

Davrinche C., Immunogenetics 32:309-312(1990).

Mejia J.E., Hum. Immunol. 39:49-53(1994).

Schwaeble W., Immunobiology 188:221-232(1993).

