



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

ITIH2 (Cleaved-Asp702) Antibody [#L0300]

Catalog Number: L0300

Concentration: 1mg/ml

Swiss-Prot No.: P19823

Other Names: Inter-alpha-trypsin inhibitor heavy chain H2; Inter-alpha-inhibitor heavy chain 2; ITI heavy chain H2; Inter-alpha-trypsin inhibitor complex component II; Serum-derived hyaluronan-associated protein; SHAP; IGHEP2

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from human ITIH2.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: ITIH2 (Cleaved-Asp702) antibody detects endogenous levels of fragment of activated ITIH2 resulting from cleavage adjacent to Asp702.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:10000

References:

Gebhard W., FEBS Lett. 229:63-67(1988).
Salier J.-P., Proc. Natl. Acad. Sci. U.S.A. 84:8272-8276(1987).
Salier J.-P., Biol. Chem. Hoppe-Seyler 369:15-18(1988).

