



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

Kallikrein-11 (Cleaved-Ile54) Antibody [#L0303]

Catalog Number: L0303

Concentration: 1mg/ml

Swiss-Prot No.: Q9UBX7

Other Names: Kallikrein-11 [Precursor]; EC 3.4.21.-; hK11; Hippostasin; Trypsin-like protease; Serine protease 20; Kallikrein-11 inactive chain 1; Kallikrein-11 inactive chain 2; KLK11; PRSS20, TLSP; UNQ649/PRO1279

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from human Kallikrein-11.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: Kallikrein-11 (Cleaved-Ile54) antibody detects endogenous levels of fragment of activated Kallikrein-11 resulting from cleavage adjacent to Ile54.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:10000

References:

Yoshida S., *Biochim. Biophys. Acta* 1399:225-228(1998).
Mitsui S., *Biochem. Biophys. Res. Commun.* 272:205-211(2000).
Yousef G.M., *Genomics* 63:88-96(2000).

