



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

MMP10 (Cleaved-Phe99) Antibody [#L0318]

Catalog Number: L0318

Concentration: 1mg/ml

Swiss-Prot No.: P09238

Other Names: STMY2; STROMELYSIN II; matrix metalloproteinase 10 (stromelysin 2)

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from human MMP10.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using

epitope-specific immunogen.

Specificity: MMP10 (Cleaved-Phe99) antibody detects endogenous levels of fragment of activated MMP10 resulting from cleavage adjacent to Phe99.

Reactivity: Human


Applications: WB: 1:500~1:3000 ELISA: 1:40000

References:

Muller D., Biochem. J. 253:187-192(1988).

Kalnine N., Submitted (OCT-2004) to the EMBL/GenBank/DDBJ databases.

Ota T., Nat. Genet. 36:40-45(2004).



For Research Use Only