



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

MMP23 (Cleaved-Tyr79) Antibody [#L0330]

Catalog Number: L0330

Concentration: 1mg/ml

Swiss-Prot No.: O75900

Other Names: Matrix metalloproteinase-23 [Precursor], EC 3.4.24.-, MMP-23, Matrix metalloproteinase 21, MMP-21, Matrix metalloproteinase 22, MMP-22, Femalysin, MIFR-1, Matrix metalloproteinase-23, soluble form, MMP23A, MMP21, MMP23B, MMP21, MMP22

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against

synthesized peptide derived from human MMP23.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: MMP23 (Cleaved-Tyr79) antibody detects endogenous levels of fragment of activated MMP23 resulting from cleavage adjacent to Tyr79.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:5000

References:

Gururajan R., Genomics 52:101-106(1998).

Velasco G., J. Biol. Chem. 274:4570-4576(1999).

Ohnishi J., Mol. Endocrinol. 15:747-764(2001).

