



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

KLKB1 (heavy chain, Cleaved-Arg390) antibody [#L0367]

Catalog Number: L0367

Concentration: 1mg/ml

Swiss-Prot No.: P03952

Other Names: EC 3.4.21.34; Fletcher factor; kininogenin; KLK3

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from human KLKB1.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: KLKB1 (heavy chain, Cleaved-Arg390) antibody detects endogenous levels of fragment of activated KLKB1 resulting from cleavage adjacent to Arg390.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:10000

References:

Chung D.W., Biochemistry 25:2410-2417(1986).

Yu H., Genomics 69:225-234(2000).

McMullen B.A., Biochemistry 30:2050-2056(1991).

