



**Web:** [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
**Order:** [order@dijibio.com](mailto:order@dijibio.com)  
**Tel:** (86) 519 8805 0026

## CHKB Antibody [#C11923]

**Catalog Number:** C11923

**Concentration:** 1mg/ml

**Swiss-Prot No.:** Q9Y259

**Other Names:** Choline/ethanolamine kinase, Choline kinase beta (EC 2.7.1.32)(CK), Ethanolamine kinase (EC 2.7.1.82)(EK), CHKB, CHETK, CHKL

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from internal of human CHKB.

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** CHKB antibody detects endogenous levels of total CHKB protein.

**Reactivity:** Human

**Applications:** WB: 1:500~1:3000 ELISA: 1:20000

**References:**

Yamazaki N., J. Biol. Chem. 275:31739-31746(2000).

Collins J.E., Genome Biol.

5:RESEARCH84.1-RESEARCH84.11(2004).

Yu L.-R., J. Proteome Res. 6:4150-4162(2007).

