



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

ASF1A Antibody [#C11928]

Catalog Number: C11928

Concentration: 1mg/ml

Swiss-Prot No.: Q9Y294

Other Names: Histone chaperone ASF1A, Anti-silencing function protein 1 homolog A, hAsf1, hAsf1a, CCG1-interacting factor A, CIA, hCIA, ASF1A, CGI-98, HSPC146

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human

ASF1A.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: ASF1A antibody detects endogenous levels of total ASF1A protein.

Reactivity: Human, Mouse

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Munakata T., Genes Cells 5:221-233(2000).

Sillje H.H.W., Curr. Biol. 11:1068-1073(2001).

Lai C.-H., Genome Res. 10:703-713(2000).

