



**Web:** [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
**Order:** [order@dijibio.com](mailto:order@dijibio.com)  
**Tel:** (86) 519 8805 0026

## ATF5 Antibody [#C11931]

**Catalog Number:** C11931

**Concentration:** 1mg/ml

**Swiss-Prot No.:** Q9Y2D1

**Other Names:** Cyclic AMP-dependent transcription factor ATF-5, Activating transcription factor 5, Transcription factor ATFx, ATF5, ATFX

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from internal of human ATF5.

**Purification:** The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** ATF5 antibody detects endogenous levels of total ATF5 protein.

**Reactivity:** Human, Mouse, Rat

**Applications:** WB: 1:500~1:3000 ELISA: 1:1000

**References:**

White J.H., Proc. Natl. Acad. Sci. U.S.A. 97:13967-13972(2000).

Yamada S., Oncogene 23:5901-5911(2004).

Liu W., Biochem. Biophys. Res. Commun. 321:954-960(2004).

