



Web: www.anbobio.com

www.dijibio.com

Order: order@dijibio.com

Tel: (86) 519 8805 0026

AKAP2 Antibody [#C11932]

Catalog Number: C11932

Concentration: 1mg/ml

Swiss-Prot No.: Q9Y2D5

Other Names: A-kinase anchor protein 2, AKA2, AKAP-2, KIAA0920, PALM2-AKAP2, PRKA2, Protein kinase A anchoring protein 2

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human

AKAP2.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: AKAP2 antibody detects endogenous levels of total AKAP2 protein.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Ota T., Nat. Genet. 36:40-45(2004).

sen J.V., Cell 127:635-648(2006).

Nagase T., DNA Res. 6:63-70(1999).

