



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

KAD6 Antibody [#C11948]

Catalog Number: C11948

Concentration: 1mg/ml

Swiss-Prot No.: Q9Y3D8

Other Names: Adenylate kinase isoenzyme 6, EC 2.7.4.3, ATP-AMP transphosphorylase 6, AK6, AD-004, CGI-137

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human KAD6.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: KAD6 antibody detects endogenous levels of total KAD6 protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Lai C.-H., Genome Res. 10:703-713(2000).

Hu R.-M., Proc. Natl. Acad. Sci. U.S.A. 97:9543-9548(2000).

Ren H., Proc. Natl. Acad. Sci. U.S.A. 102:303-308(2005).

