



Web: www.anobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

DYR1B Antibody [#C11955]

Catalog Number: C11955

Concentration: 1mg/ml

Swiss-Prot No.: Q9Y463

Other Names: DYR1B, Dual-specificity

tyrosine-phosphorylation regulated kinase 1B, EC

2.7.12.1, MIRK, Mirk protein kinase, kinase DYRK1B

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human

DYR1B.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: DYR1B antibody detects endogenous levels of total DYR1B protein.

Reactivity: Human, Mouse

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Lee K., Cancer Res. 60:3631-3637(2000).

Grimwood J., Nature 428:529-535(2004).

Lim S., . Biol. Chem. 277:25040-25046(2002).

