



**Web:** [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
**Order:** [order@dijibio.com](mailto:order@dijibio.com)  
**Tel:** (86) 519 8805 0026

## SUPT16H Antibody [#C11966]

**Catalog Number:** C11966

**Concentration:** 1mg/ml

**Swiss-Prot No.:** Q9Y5B9

**Other Names:** CDC68, Chromatin-specific transcription elongation factor 140 kDa subunit, FACT140

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from internal of human SUPT16H.

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** SUPT16H antibody detects endogenous levels of total SUPT16H protein.

**Reactivity:** Human, Mouse

**Applications:** WB: 1:500~1:3000 ELISA: 1:10000

**References:**

Kang S.-W., Genes Cells 5:251-263(2000).

Wada T., Mol. Cell 5:1067-1072(2000).

Keller D.M., Mol. Cell 7:283-292(2001).

