



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

TF3C3 Antibody [#C11971]

Catalog Number: C11971

Concentration: 1mg/ml

Swiss-Prot No.: Q9Y5Q9

Other Names: TF3C-gamma, TF3C3, TFIIIC 102 kDa subunit, TFIIIC102, TFIIICgamma, TFiiC2-102, eneral transcription factor 3C polypeptide 3, general transcription factor IIIC, polypeptide 3 (102kD), general transcription factor IIIC, polypeptide 3, 102kDa, transcription factor IIIC, 102 kD, transcription factor IIIC-gamma subunit

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against

synthesized peptide derived from internal of human TF3C3.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: TF3C3 antibody detects endogenous levels of total TF3C3 protein.

Reactivity: Human, Mouse

Applications: WB: 1:500~1:3000 ELISA: 1:1000

References:

Ota T., Nat. Genet. 36:40-45(2004).

The German cDNA consortium; Submitted (JAN-2000) to the EMBL/GenBank/DDBJ databases.

Hsieh Y.-J., Mol. Cell. Biol. 19:4944-4952(1999).

