



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

MRCKB Antibody [#C11973]

Catalog Number: C11973

Concentration: 1mg/ml

Swiss-Prot No.: Q9Y5S2

Other Names: CDC42-binding protein kinase beta, CDC42BPB, EC 2.7.11.1, MRCK-beta

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human MRCKB.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: MRCKB antibody detects endogenous levels of total MRCKB protein.

Reactivity: Human, Mouse, Rat


Applications: WB: 1:500~1:3000 IHC: 1:50~1:100
ELISA: 1:40000

References:

Nakajima D., DNA Res. 9:99-106(2002).


Olsen J.V., Cell 127:635-648(2006).

Hirosawa M., DNA Res. 6:329-336(1999).



Peptide - +

Immunohistochemistry analysis of paraffin-embedded human brain tissue using MRCKB antibody.



Western blot analysis of extracts from COLO cells, using MRCKB antibody.