



Web: www.anbio.com
www.dijbio.com
Order: order@dijbio.com
Tel: (86) 519 8805 0026

TAF6L Antibody [#C11986]

Catalog Number: C11986

Concentration: 1mg/ml

Swiss-Prot No.: Q9Y6J9

Other Names: TAF6-like RNA polymerase II
p300/CBP-associated factor-associated factor 65 kDa
subunit 6L, PCAF-associated factor 65 alpha,
PAF65-alpha, TAF6L, PAF65A

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered
saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl,
0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against
synthesized peptide derived from internal of human

TAF6L.

Purification: The antibody was affinity-purified from
rabbit antiserum by affinity-chromatography using
epitope-specific immunogen.

Specificity: TAF6L antibody detects endogenous levels
of total TAF6L protein.

Reactivity: Human, Mouse

Applications: WB: 1:500~1:3000 ELISA: 1:40000

References:

Martinez E., Mol. Cell. Biol. 21:6782-6795(2001).

Matsuoka S., Science 316:1160-1166(2007).

Struhl K., Cell 94:1-4(1998).

