



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

STON1 Antibody [#C11991]

Catalog Number: C11991

Concentration: 1mg/ml

Swiss-Prot No.: Q9Y6Q2

Other Names: SBLF, STN1; stoned

B/TFIIA-alpha/beta-like factor, STON1, stoned B-like factor, stonin-1

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from C-terminal of human

STON1.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: STON1 antibody detects endogenous levels of total STON1 protein.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:10000

References:

Martina J.A., J. Cell Biol. 153:1111-1120(2001).

Han S., Biol. Reprod. 64:507-517(2001).

Upadhyaya A.B., J. Biol. Chem. 274:18040-18048(1999).

