



Web: www.anbobio.com

www.dijibio.com

Order: order@dijibio.com

Tel: (86) 519 8805 0026

14-3-3 θ Antibody [#C12006]

Catalog Number: C12006

Concentration: 1mg/ml

Swiss-Prot No.: P27348

Other Names: 14-3-3 protein T-cell, 14-3-3 protein tau, 1433T, 143T, HS1 protein, YWHAQ

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human 14-3-3 θ .

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: 14-3-3 θ antibody detects endogenous levels of total 14-3-3 θ protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Nielsen P.J., Biochim. Biophys. Acta 1088:425-428(1991).

Leffers H., J. Mol. Biol. 231:982-998(1993).

Gevaert K., Nat. Biotechnol. 21:566-569(2003).

