



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

5-HT-1F Antibody [#C12012]

Catalog Number: C12012

Concentration: 1mg/ml

Swiss-Prot No.: P30939

Other Names: 5-hydroxytryptamine receptor 1F, 5-HT-1F, Serotonin receptor 1F, HTR1F, HTR1EL

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human 5-HT-1F.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: 5-HT-1F antibody detects endogenous levels of total 5-HT-1F protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:40000

References:

Lovenberg T.W., Proc. Natl. Acad. Sci. U.S.A. 90:2184-2188(1993).

Adham N., Proc. Natl. Acad. Sci. U.S.A. 90:408-412(1993).

The MGC Project Team, Genome Res. 14:2121-2127(2004).

