



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

5-HT-3A Antibody [#C12016]

Catalog Number: C12016

Concentration: 1mg/ml

Swiss-Prot No.: P46098

Other Names: 5-HT-3, 5-HT3R, 5-hydroxytryptamine 3 receptor, 5HT3R, HTR3, HTR3A, serotonin-gated ion channel receptor

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human

5-HT-3A.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: 5-HT-3A antibody detects endogenous levels of total 5-HT-3A protein.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:1000

References:

Miyake A., Mol. Pharmacol. 48:407-416(1995).

Belelli D., Mol. Pharmacol. 48:1054-1062(1995).

Yamada K., Biol. Psychiatry 60:192-201(2006).

