



**Web:** [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
**Order:** [order@dijibio.com](mailto:order@dijibio.com)  
**Tel:** (86) 519 8805 0026

## 5-HT-6 Antibody [#C12019]

**Catalog Number:** C12019

**Concentration:** 1mg/ml

**Swiss-Prot No.:** P50406

**Other Names:** 5-HT6; 5HT6R; HTR6

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from internal of human 5-HT-6.

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using

epitope-specific immunogen.

**Specificity:** 5-HT-6 antibody detects endogenous levels of total 5-HT-6 protein.

**Reactivity:** Human (Identities = 100%, Positives = 100%);  
Mouse (Identities = 100%, Positives = 100%);  
Rat (Identities = 100%, Positives = 100%)

**Applications:** WB: 1:500~1:3000 ELISA: 1:20000

**References:**

Kohen R., J. Neurochem. 66:47-56(1996).

Kopatz S.A., Submitted (OCT-2003) to the EMBL/GenBank/DDBJ databases.

Gregory S.G., Nature 441:315-321(2006).

