



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

ACTN4 Antibody [#C12027]

Catalog Number: C12027

Concentration: 1mg/ml

Swiss-Prot No.: P12814/P35609/Q08043/O43707

Other Names: AAC4, Alpha-actinin 4, F-actin cross linking protein, Non-muscle alpha-actinin 4, actinin, alpha 4

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human

ACTN4.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: ACTN4 antibody detects endogenous levels of total ACTN4 protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Honda K., J. Cell Biol. 140:1383-1393(1998).

Nikolopoulos S.N., gene 19:380-386(2000).

Patrie K.M., . Biol. Chem. 277:30183-30190(2002).

