



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

ADCY1 Antibody [#C12031]

Catalog Number: C12031

Concentration: 1mg/ml

Swiss-Prot No.: Q08828

Other Names: ATP pyrophosphate-lyase, CYA1, Ca/calmodulin activated adenylyl cyclase, EC 4.6.1.1, adenylate cyclase 1, adenylate cyclase 1 (brain), adenylate cyclase, type I

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against

synthesized peptide derived from internal of human ADCY1.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: ADCY1 antibody detects endogenous levels of total ADCY1 protein.

Reactivity: Human, Mouse

Applications: WB: 1:500~1:3000 ELISA: 1:40000

References:

Hillier L.W., Nature 424:157-164(2003).

Villacres E.C., Genomics 16:473-478(1993).

