



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

ADRB2 Antibody [#C12044]

Catalog Number: C12044

Concentration: 1mg/ml

Swiss-Prot No.: P07550

Other Names: ADRB2R, B2AR, adrenergic receptor, beta-2, beta-2 adrenergic receptor

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from C-terminal of human ADRB2.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: ADRB2 antibody detects endogenous levels of total ADRB2 protein.

Reactivity: Human, Mouse, Rat


Applications: WB: 1:500~1:3000 IF: 1:100~1:500
ELISA: 1:20000

References:


Moffett S., J. Neurochem. 76:269-279(2001).

Cao T.T., Nature 401:286-290(1999).

Whistler J.L., Science 297:615-620(2002).



Western blot analysis of extracts from HeLa cells, using ADRB2 antibody.



Immunofluorescence analysis of NIH/3T3 cells, using ADRB2 antibody.

For Research Use Only