



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

AIFM3 Antibody [#C12048]

Catalog Number: C12048

Concentration: 1mg/ml

Swiss-Prot No.: Q96NN9

Other Names: Apoptosis-inducing factor 3, EC 1.-.-.-, Apoptosis-inducing factor-like protein, AIFM3, AIFL

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human AIFM3.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: AIFM3 antibody detects endogenous levels of total AIFM3 protein.

Reactivity: Human (Identities = 100%, Positives = 100%);
Mouse (Identities = 100%, Positives = 100%)


Applications: IHC: 1:50~1:100 IF: 1:100~1:500
ELISA: 1:40000

References:

Xie Q., J. Biol. Chem. 280:19673-19681(2005).


Collins J.E., Genome Biol. 5:RESEARCH84.1-RESEARCH84.11(2004).

Ota T., Nat. Genet. 36:40-45(2004).



Peptide - +

Immunohistochemistry analysis of paraffin-embedded human brain tissue, using AIFM3 antibody.



Peptide - +

Immunofluorescence analysis of A549 cells, using AIFM3 antibody.