



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

API-5 Antibody [#C12052]

Catalog Number: C12052

Concentration: 1mg/ml

Swiss-Prot No.: Q9BZZ5

Other Names: AAC-11, API-5, API5-like 1, API5L1, FIF, XAGL protein, apoptosis inhibitor 5, fibroblast growth factor 2-interacting factor, fibroblast growth factor 2-interacting factor 2, migration-inducing protein MIG8

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from C-terminal of human

API-5.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: API-5 antibody detects endogenous levels of total API-5 protein.

Reactivity: Human, Mouse

Applications: WB: 1:500~1:3000 IHC: 1:50~1:100


ELISA: 1:20000

References:

Tewari M., Cancer Res. 57:4063-4069(1997).

Van den Berghe L., Mol. Endocrinol. 14:1709-1724(2000).


Olsen J.V., Cell 127:635-648(2006).



Peptide -

+

Immunohistochemistry analysis of paraffin-embedded human heart tissue using API-5 antibody.



Western blot analysis of extracts from 293 cells, Jurkat cells and HUVEC cells, using API-5 antibody.