



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

ARRDC3 Antibody [#C12061]

Catalog Number: C12061

Concentration: 1mg/ml

Swiss-Prot No.: Q96B67

Other Names: Arrestin domain-containing protein 3, ARRDC3, KIAA1376

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from C-terminal of human ARRDC3.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: ARRDC3 antibody detects endogenous levels of total ARRDC3 protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Nagase T., DNA Res. 7:65-73(2000).

The MGC Project Team, Genome Res. 14:2121-2127(2004).

