



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

BUB3 Antibody [#C12075]

Catalog Number: C12075

Concentration: 1mg/ml

Swiss-Prot No.: O43684

Other Names: Mitotic checkpoint protein BUB3, BUB3

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from C-terminal of human BUB3.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: BUB3 antibody detects endogenous levels of total BUB3 protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:1000

References:

Taylor S.S., J. Cell Biol. 142:1-11(1998).

The MGC Project Team, Genome Res. 14:2121-2127(2004).

Beausoleil S.A., Proc. Natl. Acad. Sci. U.S.A. 101:12130-12135(2004).

