



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

C56D2 Antibody [#C12082]

Catalog Number: C12082

Concentration: 1mg/ml

Swiss-Prot No.: O14569

Other Names: Cytochrome b561 domain-containing protein 2; Putative tumor suppressor protein 101F6; CYB561D2; 101F6; LUCA12.2

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human C56D2.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: C56D2 antibody detects endogenous levels of total C56D2 protein.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Bader S., Submitted (JAN-1998) to the EMBL/GenBank/DDBJ databases.

Muzny D.M., Nature 440:1194-1198(2006).

The MGC Project Team; Genome Res. 14:2121-2127(2004).

