



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

CDH15 Antibody [#C12090]

Catalog Number: C12090

Concentration: 1mg/ml

Swiss-Prot No.: P55291

Other Names: Cadherin-15 [Precursor], Muscle cadherin, M-cadherin, Cadherin-14, CDH15, CDH14, CDH3

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human CDH15.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: CDH15 antibody detects endogenous levels of total CDH15 protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:10000

References:

Shibata T., J. Biol. Chem. 272:5236-5240(1997).

The MGC Project Team, Genome Res. 14:2121-2127(2004).

