



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

CDH16 Antibody [#C12091]

Catalog Number: C12091

Concentration: 1mg/ml

Swiss-Prot No.: O75309

Other Names: Cadherin-16 [Precursor], Kidney-specific cadherin, Ksp-cadherin, CDH16, UNQ695/PRO1340

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from C-terminal of human CDH16.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: CDH16 antibody detects endogenous levels of total CDH16 protein.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:10000

References:

Thomson R.B., Genomics 51:445-451(1998).

Clark H.F., Genome Res. 13:2265-2270(2003).

The MGC Project Team, Genome Res. 14:2121-2127(2004).

