



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

CDH18 Antibody [#C12093]

Catalog Number: C12093

Concentration: 1mg/ml

Swiss-Prot No.: Q13634

Other Names: Cadherin-18 [Precursor], Cadherin-14, CDH18, CDH14

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human

CDH18.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: CDH18 antibody detects endogenous levels of total CDH18 protein.

Reactivity: Human, Mouse

Applications: WB: 1:500~1:3000 ELISA: 1:10000

References:

Shibata T., J. Biol. Chem. 272:5236-5240(1997).

