



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

CDH19 Antibody [#C12094]

Catalog Number: C12094

Concentration: 1mg/ml

Swiss-Prot No.: Q9H159

Other Names: Cadherin-19 [Precursor], CDH19, CDH7L2, UNQ478/PRO941, CAD19

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human

CDH19.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: CDH19 antibody detects endogenous levels of total CDH19 protein.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Kools P., Genomics 68:283-295(2000).

Clark H.F., Genome Res. 13:2265-2270(2003).

