



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

CDH2 Antibody [#C12100]

Catalog Number: C12100

Concentration: 1mg/ml

Swiss-Prot No.: P19022

Other Names: CADH2, CDHN, Cadherin-2, N-cad, N-cadherin, NCAD, Neural-cadherin precursor, cadherin, neural

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human

CDH2.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: CDH2 antibody detects endogenous levels of total CDH2 protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:10000

References:

Reid R.A., Nucleic Acids Res. 18:5896-5896(1990).

Salomon D., J. Cell Sci. 102:7-17(1992).

Amanchy R., J. Proteome Res. 4:1661-1671(2005).

