



**Web:** [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
**Order:** [order@dijibio.com](mailto:order@dijibio.com)  
**Tel:** (86) 519 8805 0026

## TMEM30B Antibody [#C12119]

**Catalog Number:** C12119

**Concentration:** 1mg/ml

**Swiss-Prot No.:** Q3MIR4

**Other Names:** Cell cycle control protein 50B,  
Transmembrane protein 30B, TMEM30B, CDC50B,  
CC50B

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered  
saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl,  
0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against  
synthesized peptide derived from internal of human

TMEM30B.

**Purification:** The antibody was affinity-purified from  
rabbit antiserum by affinity-chromatography using  
epitope-specific immunogen.

**Specificity:** TMEM30B antibody detects endogenous  
levels of total TMEM30B protein.

**Reactivity:** Human (Identities = 100%, Positives = 100%)

**Applications:** WB: 1:500~1:3000 ELISA: 1:40000

**References:**

The MGC Project Team, Genome Res.  
14:2121-2127(2004).

Katoh Y., Oncol. Rep. 12:939-943(2004).

