



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

CDCA4 Antibody [#C12175]

Catalog Number: C12175

Concentration: 1mg/ml

Swiss-Prot No.: Q9BXL8

Other Names: Cell division cycle-associated protein 4, Hematopoietic progenitor protein, HEPP

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human CDCA4.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: CDCA4 antibody detects endogenous levels of total CDCA4 protein.


Reactivity: Human (Identities = 100%, Positives = 100%);
Mouse (Identities = 83%, Positives = 91%)

Applications: IHC: 1:50~1:100 IF: 1:100~1:500
ELISA: 1:1000

References:


Abdullah J.M., Blood Cells Mol. Dis. 27:667-676(2001).

Ota T., Suzuki Y., Nishikawa T., Nat. Genet. 36:40-45(2004).



Peptide - +

Immunohistochemistry analysis of paraffin-embedded human brain tissue using CDCA4 antibody.



Peptide - +

Immunofluorescence analysis of HepG2 cells, using CDCA4 antibody.