



**Web:** [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
**Order:** [order@dijibio.com](mailto:order@dijibio.com)  
**Tel:** (86) 519 8805 0026

## CDKAL Antibody [#C12178]

**Catalog Number:** C12178

**Concentration:** 1mg/ml

**Swiss-Prot No.:** Q5VV42

**Other Names:** CDK5 regulatory subunit-associated protein 1-like 1; CDKAL

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from N-terminal of human CDKAL.

**Purification:** The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** CDKAL antibody detects endogenous levels of total CDKAL protein.

**Reactivity:** Human (Identities = 100%, Positives = 100%);  
Rat (Identities = 92%, Positives = 100%)

**Applications:** WB: 1:500~1:3000 ELISA: 1:40000

**References:**

Ota T., Nat. Genet. 36:40-45(2004).

Mungall A.J., Nature 425:805-811(2003).

Mural R.J., Submitted (JUL-2005) to the EMBL/GenBank/DDBJ databases.

