



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

RXFP3 Antibody [#G738]

Catalog Number: G738

Concentration: 1mg/ml

Swiss-Prot No.: Q9NSD7

Other Names: G-protein coupled receptor SALPR; GPCR135; relaxin-3 receptor 1; RL3R1; SALPR; somatostatin- and angiotensin-like peptide receptor

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human RXFP3.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: RXFP3 antibody detects endogenous levels of total RXFP3 protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:10000

References:

Matsumoto M., Gene 248:183-189(2000).

Kopatz S.A., Submitted (FEB-2003) to the EMBL/GenBank/DDBJ databases.

Liu C., J. Biol. Chem. 278:50754-50764(2003).

