



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

SCTR Antibody [#G742]

Catalog Number: G742

Concentration: 1mg/ml

Swiss-Prot No.: P47872

Other Names: Secretin receptor [Precursor]; SCT-R

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human SCTR.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: SCTR antibody detects endogenous levels of total SCTR protein.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Jiang S., Biochem. Biophys. Res. Commun. 207:883-890(1995).

Chow B.K.-C.; Biochem. Biophys. Res. Commun. 212:204-211(1995).

Patel D.R., Mol. Pharmacol. 47:467-473(1995).

