



**Web:** [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
**Order:** [order@dijibio.com](mailto:order@dijibio.com)  
**Tel:** (86) 519 8805 0026

## SSTR4 Antibody [#G747]

**Catalog Number:** G747

**Concentration:** 1mg/ml

**Swiss-Prot No.:** P31391

**Other Names:** G-protein coupled receptor; somatostatin receptor 4; SS4R; SSR4

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from internal of human SSTR4.

**Purification:** The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** SSTR4 antibody detects endogenous levels of total SSTR4 protein.

**Reactivity:** Human, Mouse, Rat

**Applications:** WB: 1:500~1:3000 ELISA: 1:20000

**References:**

Xu Y., Biochem. Biophys. Res. Commun. 193:648-652(1993).

Rohrer L., c. Natl. Acad. Sci. U.S.A. 90:4196-4200(1993).

Yamada Y., Biochem. Biophys. Res. Commun. 195:844-852(1993).

